

Sperm search after using SpermElute™ is five times more sensitive than the P30 assay

Sample	Dilution	Target Sperm added	Method	P30	Sperm Search Positive
396	1 to 1000	465	VORTEX	POS	POS (2)
397	1 to 2500	251	VORTEX	NEG	NEG
398	1 to 5000	137	VORTEX	NEG	NEG
399	1 to 7500	43	VORTEX	NEG	NEG
400	1 to 10000	21	VORTEX	NEG	NEG
401	1 to 1000	465	MO-Lite	POS	POS (15)
402	1 to 2500	251	MO-Lite	NEG	POS (8)
403	1 to 5000	137	MO-Lite	NEG	POS (4)
404	1 to 7500	43	MO-Lite	NEG	NEG
405	1 to 10000	21	MO-Lite	NEG	NEG